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TEMPORARY TRANSVENOUS PACING ELECTRODE

Type: ES and ESM

Product compliance with Directive 93/42/EEC, number of the Notified Body


HAGMED Sp. z o.o. sp. k.
96-200 Rawa Mazowiecka ul. Tomaszowska 32; POLAND
Tel. +48 46 814 44 29; Fax. +48 46 814 48 05
e-mail info@hagmed.com www.hagmed.com

INSTRUCTION FOR USE

1. PRODUCT DESCRIPTION AND INDICATION FOR USE.

The temporary transvenous pacing electrode is sterile, single use medical device. The electrode is used for temporary transvenous cardiac pacing (up to 2 weeks) mainly for reanimation of patients with sudden ventricular arrest. It can be used in current or threatening bradycardia caused by conduction system disease or acute myocardial infarction. Indication for use may be treatments, during which may occur clinically significant decrease of heart rate, as well as the need to control recurrent, dangerous ventricular arrhythmias.


Basic technical specification and catalogue numbers "REF" was placed in the product catalogue. Information identifying the product in accordance with catalogue data are placed on the product label. Essential parts of the electrode and its accessories are shown in the below figure.

1 – Distal ring (-), 2 – Proximal ring (+), 3 – Main shaft, 4 – Central connection, 5 – PROX plug (+), 6 – DISTAL plug (-), 7 – PROX plug shroud, 8 – DISTAL plug shroud

Parts 1, 2 and 3 form the electrode working part. The electrode working part is radiologically detectable.

The electrode is scaled at 10 cm according to the below figure:

2. CONTRAINDICATIONS.

Lack of the patient's consent is absolute contraindication to use the temporary transvenous pacing electrode. Relative contraindications are situations amenable to modification or possible to ignore, in a case of expected greater benefits against risks arising from relative contraindications. To relative contraindications should be included: implanted artificial tricuspid valve, atria tricuspid valve, high risk of central venous cannulation, blood coagulation disorders, active or general recurrent infections, unstable patients exposed to occurrence of serious arrhythmias caused by introduction of the electrode. Other contraindications to the temporary transvenous cardiac pacing may also occur, the final decision to commence the procedure belongs to the physician.

3. RECOMMENDATIONS FOR USE AND PRECAUTIONS.

The temporary transvenous pacing electrode is recommended to be inserted into the right ventricle or the right atrium of the heart from the cervical, subclavian or femoral access by puncturing a blood vessel with use vascular cannulas or introducer kits, according to instruction for use of used vascular access equipment.

PREPARING OF ELECTRODE TO INSERTION PROCEDURE

Check expiration date and potential damage of the electrode single packaging. Do not use the electrode if expiration date has been exceeded or the single packaging is damaged.

Remove the electrode from the single packaging and place in the sterile operation field, directly before use.

Check carefully technical condition of the electrode. Do not use of the electrode if the distal end or the main shaft were deformed, fragmented or damaged in other perceptible way. Avoid manoeuvres that may damage the electrode, catch with use sharp-edge tools, stretch, coil, etc. Use of organic solvents on the surfaces of electrode elements is forbidden. Do not allow moisture the electrode plugs.

Check dimensional compatibility of the electrode and the introducer kit or the vascular cannula.

Selection of vascular cannulas or introducers kits depending on the electrode size			List of external pacemakers, compatible with the electrodes Biotronic: EDP 20/A; EDP 20/B; EDP 30/A; EDP 30/B; Recor S; Recor D St. Jude Medical: 3085; 5391 Medtronic: 5391; 5348; 5388 ITAM: MIP-801
Electrode size	Cannula size	Introducer size	
4F	G16 or G14	4F or 5F	
5F	G14 or G12	5F or 6F	
6F	G12	6F or 7F	
7F	Is not recommend	7F or 8F	

INSERTION PROCEDURE OF ELECTRODE


Before proceeding to insert the electrode into the cardiovascular system, it is recommended to check the electrode's electrical installations with use appropriate control functions of used external pacemaker, in accordance with its instruction for use. In the event of failure, the electrode shall be replaced.

Carefully and slowly insert the electrode inside the circulatory system and the heart with use the electrocardiographic (ECG) monitoring. Avoid excessive force during inserting or withdrawing of the electrode.

During inserting the electrode with use ECG monitoring, the DISTAL plug is connected to defined lead (V1, I or II) of the electrocardiograph, according to its instruction for use.

During inserting of the electrode, ECG records should be observed in the following sequence: peripheral vessels, the superior vena cava, the right atrium and the right ventricle. It is recommended to verify the positioning of the electrode distal end in the heart with use the X-ray imaging.

After obtaining a satisfactory location of the electrode distal end in the heart, remove the cannula or the introducer from a blood vessel. If cannulas or introducers "peel away" type are used, the process of they withdrawal should be gradual and connected with simultaneous tearing with use the "T" connectors. Tearing process must not take place in a blood vessel.

Attach the electrode to the patient's skin in order to prevent dislocation and then perform a sterile dressing of the puncture site.

Determine the stimulation threshold and then connect the electrode to the external pacemaker in the following way: the DISTAL electrode plug with the external pacemaker negative socket, the PROX electrode plug with the external pacemaker positive socket. The connection of the electrode with the external pacemaker can also be performed with use an extension cable being an equipment of the external pacemaker.

Stimulation parameters should be set depending on the patient clinical situation.

Process of the electrode withdrawal should be carried out with use the ECG monitoring. After removal of the electrode from the patient's body, the electrode technical condition shall be checked in order to exclude fragmentation or potential lodged of the electrode fragments in the patient's heart.

For recording electrophysiological parameters from the electrode and conducting temporary heart pacing, only compatible electrocardiographs and external pacemakers with CF protection type (marked with the symbol "heart in square"), having insulated patient circuit, shall be used.

Leakage current from the electrocardiograph to which the electrode has been connected must not exceed 10µA.

Do not leave unprotected electrode plugs, to avoid uncontrolled flow to the patient's heart of electrical charge from others devices, causing ventricular fibrillation.

During preparation and realization of the temporary transvenous cardiac pacing, medical procedures and complete asepsis shall be maintained.

The temporary transvenous pacing electrode shall only be used by medical personnel with appropriate qualifications and experience in the electrophysiology field and use of intracardiac electrodes.

The temporary transvenous pacing electrode after use shall be subjected to a recycling process. Re-sterilization and re-use of the temporary transvenous pacing electrode is forbidden.

4. POSSIBLE RISKS AND COMPLICATIONS.

Following risks and complications are associated with use of the temporary transvenous pacing electrode: vascular access infection, hematoma at vascular access site, access site vessel occlusion, arteriovenous fistula, deep vein thrombosis, air embolism, vessel perforation, pleural hematoma, tricuspid valve damage, cardiac wall perforation, atrio-ventricular fistula, cardiac tamponade, pericardial effusion, lung perforation, hypotension, pneumothorax, hematoma in heart cavities, electrode trapping in cardiac structures, electrode dislocation, electrode fragmentation in the heart, transient or permanent atrioventricular block, unintentional induction of unwanted cardiac arrhythmia, including asystole, thromboembolic complications, including pulmonary embolism and brain stroke, sepsis, death.

5. TRANSPORT AND STORAGE.

The permissible storage temperature of the temporary transvenous pacing electrode is from +5°C to +45°C. The permissible transport temperature of the electrode is from -10°C to +50°C. The electrode should be protected from moisture and sunlight.




6. RESPONSIBILITY.

HAGMED Sp. z o.o. sp. k. is not liable for damages resulting from use of the product contrary to its purpose or its use contrary to recommendations of the product instruction for use. The manufacturer is not liable for damages resulting from use of the re-sterilized product or modified on their own, or use of the product with damaged barrier packing, e.g. during transport, storage, etc. Taking into account biological differences of patients, the product effectiveness cannot be completely guaranteed.

If during use, a defect of the product has been detected, the defective product must be preserved and protected with its original packaging. HAGMED Sp. z o.o. sp. k. must be informed about the incident.

The original version of this instruction for use was developed in the Polish language. In case of a problem with interpretation of translated versions, please refer to text of the original version.

Edition: ES/1/2020/EN

Description of labels symbols		
REF	Catalogue number	 Manufacturer
LOT	LOT number	 Type CF applied part
STERILE R	Sterilized using irradiation	 Non-pyrogenic
	Use by date	 Do not use if package is damaged
	Consult instruction for use	 Keep away from sunlight
	Do not re-use	 Keep dry
	Do not re-sterilize	 Temperature limit